


MultiLab® Series 2

Series 2-LHS TI

Unetixs puts the Volume back in Pulse Volume Recording. Calibrated PVR allows definitive post intervention follow-up with diagnostic quality waveforms.

Accuracy is guaranteed with the Unetixs Doppler, unsurpassed in signal strength and tonal quality.

Our easily customized procedures satisfy both the highly trained vascular technologist and the new user.



Series 2-LHS IMG

Save financial and facility resources with the only combination physiologic and fully-featured duplex system available.

Full Reporting Package

Superior Color Resolution with Triplex

Unetixs Signature Touch Screen Controls

MultiLab® Series 2-CP

A fully capable dual channel system with built in CPU and mass storage in a truly portable package.

Perfect mid-range system for those who will not compromise on quality results.

Touch Screen Controls, Full-Function Electronic Manifold and Automatic Inflation will please technologists of all skill levels.



A Truly Portable Unit

Optional Rolling Carrying Case

Optional 2-CP System Cart

The 2-CP is similar to the LHS TI in all respects except for size and level of automation.

Revo® 1100

The solution for diagnosing one of the biggest health problems in the world today is remarkably small.

Revo® is the One Touch PAD Solution.



Featuring TruDop®

Doppler Probes

and the Unetixs PPG Digit Clips and Sensors

- Automated Diagnostic ABI in as little as 5 Minutes
- Direct Printing of Full Page Color Reports
- Dual Channel C-PVR™ (Calibrated Pulse Volume Recording)
- Dual Channel PPG (Photoplethysmography) and Pressures
- Composite TruDop® CW Doppler Probes
- Automated 4-Blood Pressure Cuff Manifold
- Full Color, High Resolution Touch Screen
- Windows XP Operating System

Patent Pending