

MULTILAB® SERIES II  
LHS

GENERAL SPECIFICATIONS

Instrument Features:

- 17" Flat Touch Screen Color Monitor
- Pentium or Equivalent CPU w/ Hard Disk Drive
- DVD Writer and 3½" Floppy Drive
- Windows XP Operating System
- Color Inkjet Printer (Hewlett Packard usually supplied, HP Laser Printer optional)
- Custom Design Cart
- Built-in Hospital Isolation Module with Main Power Switch
- Stereo Headphones (customer option)
- 14 Unetixs Pressure Cuffs (assorted sizes)
- Stereo Doppler
- Full-function Remote Control
- Bilateral Testing Capabilities
- Custom Software-configurable Test Protocols and Reports

The MultiLab® Series II LHS is designed to perform all standard noninvasive arterial and venous studies incorporating:

- 8 MHz Bidirectional Doppler Probe
- 5 MHz Bidirectional Doppler Probe
- Calibrated Dual Channel Pulse Volume Recording (C-PVR™)
- Dual Channel Photoplethysmography

Electronic Cuff Manifold Features:

- 10-Cuff Simultaneous Hookup
- Full-function Remote Control for Automatic Inflation, Deflation, Capture, Forward, Back, Scrolling and more

Cart Dimensions:

- Base Size 20" w x 24" d x 40" h (does not include monitor height)
- Latching Storage Drawer, 18" w x 6" h x 20" d
- Printer Shelf
- Tower CPU Compartment
- 4-way Locking Casters with True Track
- Approx. Weight Complete: 155lbs.
- Power: AC only, 115V 50/60Hz (Note: can be supplied with other voltages)
- Integrated Full-function Keyboard

Design and system specifications subject to change/enhancement without notice. System components and accessories may not appear as shown.

The following are trademarks of their respective owners:  
Deskjet of Hewlett Packard Corporation, Microsoft Windows XP



Interpretation Manual for Basic Physiologic Testing, MultiLab Series II Operator's Manual and items pictured above supplied in the Standard Configuration

★ Infrared Temperature Instrument optional on both LHS and LHS-TI models

WWW.UNETIXS.COM  
(800)486-3849 sales@unetixs.com  
Tel: (401)294-7559 Fax: (401)294-3893  
International inquiries contact: international@unetixs.com

MULTILAB® SERIES II  
LHS

NEW  
FEATURES FOR  
MULTILAB® SYSTEMS!  
(SEE BELOW\*)



\*NOW FEATURING

TruDop®




EXCLUSIVE  
COMPOSITE TIP  
TECHNOLOGY  
5 OR 8MHz

THE NEW STANDARD FOR  
CW DOPPLER PROBE ACCURACY, DURABILITY  
AND COMFORT Compatible with Unetixs MultiLab® Systems only.

AND THE NEW UNETIXS

PPG DIGIT  
CLIPS AND  
SENSORS



Features and Benefits:

- Speed - No application of tape or Velcro needed for most studies utilizing the PPG.
- Uniform Application - Same pressure applied bilaterally.
- Sanitary - The sensors and clip can be easily cleaned.
- Flexible - Accommodates all patient types, from the very large to the very small

Design and system specifications subject to change/enhancement without notice. System components and accessories may not appear as shown.

WWW.UNETIXS.COM  
(800)486-3849 sales@unetixs.com  
International inquiries contact: international@unetixs.com

