

VASCULINK Worksheet



A. Information from the Conformance Statement for the client's SCP:

- 1. Ensure Secondary Capture Storage SOP Class is supported (1.2.840.10008.5.1.4.1.1)? YES NO
- 2. is Study identified by Accession number? YES NO
- 3. If answered "No" to above, is Study identified by other? Specify: _____

B. Information from the client's IT Department or DICOM Administrator:

- 1. MultiLab Revo IP address: _____ DHCP
- 2. MultiLab Revo Gateway: _____ DHCP
- 3. MultiLab Revo Subnet Mask: _____ DHCP
- 4. MultiLab Revo AE Title: _____
- 5. Store SCP (PACS server) AE Title: _____
- 6. Store SCP (PACS server) IP address: _____
- 7. Store SCP port number (port the SCP receives on): _____
- 8. Worklist SCP AE Title: (Worklist not currently supported by REVO) _____
- 9. Worklist SCP IP address: (Worklist not currently supported by REVO) _____
- 10. Worklist SCP port number: (Worklist not currently supported by REVO) _____
- 11. PACS Vendor and Software Version: _____